



PERSON-PORTABLE GC-MS CHEMICAL IDENTIFIER

FLIR GRIFFIN™ G510

The FLIR Griffin™ G510 Gas Chromatograph Mass Spectrometer (GC/MS) is a versatile, person-portable chemical identifier. It complements presumptive techniques used during emergency missions, by enabling responders to analyze all phases of matter (liquid, solid, vapor) and by performing rapid field- confirmation of chemical hazards. The integrated heated sample probe enables hot zone operators to identify vapor-phase chemical threats within seconds when operated in Survey Mode. The integrated split/splitless injector allows for environmental, forensic, and hazardous material sampling via syringe injection of organic liquids. The 9″ on-board touchscreen displays guided user prompts and can be operated while wearing full personal protective equipment (PPE) downrange. It is built with an IP65-rated enclosure for harsh environments and supports passive defense, interdiction, elimination, and consequence management missions. Long-lasting, on-board batteries ensure every mission is supported from beginning to end.

www.flir.com/G510



MASS SPEC PERFORMANCE REDEFINED

Confidently identify unknowns and take action with guided prompts and simple threat alarms

- Laboratory level, gold-standard linear quadrupole mass analyzer
- Full NIST and SWGDRUG on-board libraries provide confirmatory identification and analysis of trace compounds, unknown chemicals, and mixtures
- Simple on-board touchscreen with navigation assistant and Method Selector tool
- Visual and audible alarm confirmation with limited data interpretation
- On-board WiFi and GPS assist in maintaining result defensibility
- Remote monitoring and control with WiFi hotspot



ULTIMATE CHEMICAL DETECTION TOOLBOX

Versatile in-field sampling options for vapor, liquid, and solid samples

- Vapor sampling probe with rapid-response
 survey mode
- Integrated split/splitless liquid injector accepts direct injection of organic liquids
- Available Prepless Sample Introduction (PSI) Probe, Sample Prep Kit (SPK), and Touch-And-Go (TAG™) capability for direct analysis of solid samples
- Effortlessly links with SPME and headspace sample collection tools
- High-fidelity, low thermal mass (LTM) GC column for unsurpassed resolution in challenging environments



BUILT FOR EVERYONE, BUILT FOR EVERYWHERE

Completely self-contained and mission-ready from the field to the lab

- IP65-rated, dust-tight and spray-resistant
- Built-in active pumping system eliminates need for an external service module
- Integrated carrier gas, batteries, and training reference videos
- Simple field maintenance activities for increased uptime
- Extensive training, service, and reachback options available
- Optional vehicle mount kit for shock and vibration protection during off-road or harsh transport

SPECIFICATIONS

Griffin G510

Technology

Weight

Dimensions (L x W x H)

Gas Chromatography/Mass Spectrometry (GC/MS)

13.25 x 13.25 x 15.75 in (33.7 x 33.7 x 40 cm) - includes batteries, carrier gas, and vacuum system

36 lbs (16.3 kg) - includes batteries, carrier gas, and vacuum system

32 to 104 °F (0 to 40 °C); <95% relative humidity

-13 to 131 °F (-25 to 55 °C)

Road-Ready Version: 16.7 x 14 x 15.2 in (with G510 installed)

Optional Vehicle Mount Kit (LxWxH)

Operating Temp / Humidity

Storage Temp

Decontamination

Power Supply

Battery Life

Start Up Time

Carrier Gas

Calibrant

System Interface

Display Alerts

9" Multitouch Color Display

(1280x720 WVGA;1300 nits brightness)

Audible and visual (touchscreen and handheld probe)

GSS Level 1 Touch; multiple user levels Software

Communication 2 x USB 2.0, Bluetooth 4.0, WiFi 802.11n, Ethernet via USB,

integrated GPS

Internal ≥128GB SSD Data Storage

Training Requirements 2 hours basic operation; 8 hours Operator Certification

Sampling & Identification

Sample Phase

Threats

Sample Introduction

Solid, liquid, and vapor

Heated Sample Probe (included standard):

- Vapor survey mode via Membrane Introduction Mass Spectrometry (MIMS) Inlet
- Vapor confirmation via Internal Dual-Bed Preconcentrator

Split/splitless injector (included standard) accepts:

- Direct liquid sampling (organic solution) via syringe
- Liquid extraction via SPME fiber or PSI-Probe w/ Gerstel Twister™
- Solid PSI-Probe™ thermal separation via TAG™ *

Sample Prep Kit (SPK)

Simplified analysis of solids and liquids

Detects and identifies trace-level explosives, narcotics like Fentanyl and its analogues, A-series and Novichok CWAs, TICs,

environmental pollutants, and other chemicals

Standard Reference Database NIST/EPA/NIH Mass Spectral Library, SWGDRUG Mass Spectral Library, and GriffinLib Mass Spectral Library included,

with instant updates available

*optional accessories



Full identification in 4-15 minutes for most chemicals; identification within seconds (near real-time) when operating in Survey Mode

Electron Impact Ionization; non-radioactive ionization source

Self-contained miniature turbomolecular & diaphragm pumps

PPM (parts per million) - PPT (parts per trillion)

DB-5MS (15 m x .18 mm x .18 um); others available

Programmable 40 to 300 °C; ramping of 100 °C/min

Linear quadrupole mass filter

15-515 m/z; 0.7 amu @ FWHM

Electron Multiplier

7 decades

Mass Spectrometer

Mass Analyzer Type

Mass Range / Resolution

Ionization Type / Source

Vacuum System

Dynamic Range

LTM-GC Column

Temperature Range

Detection Limit Gas Chromatograph 100-240V 50-60Hz (220 W max); 19V (DC); 2 x #2590 @ 15V Li

Ion batteries (included)

Sealed for Survey Mode operation in hot-zone; IP65-rated

Off-Road Version: 16.7 x 17.5 x 16.2 in (with G510 installed)

4 hrs in Survey Mode, 2 hrs in Confirmation Mode;

hot swappable 15 minutes to full operation from cold

enclosure is dust-tight and spray-resistant

On-board FC-43 (Perfluorotributylamine)

On-board helium: external helium connector, automatic switching (Hydrogen capable)

Specifications are subject to change without notice. For the most up-to-date specs, go to www.flir.com





Method Selector Tool



Trace-level detection of narcotics like Fentanyl and its analogues



GSS Touch Chemical Identification Result

HEADQUARTERS

FLIR Systems, Inc. 27700 SW Parkway Ave. Wilsonville, OR 97070 USA

DETECTION SALES, AMERICAS

FLIR Systems, Inc. 1201 S. Joyce Street Suite C006 Arlington, VA 22202 PH: +1-877-692-2120

DETECTION SALES, APAC

FLIR Detection, Inc. 10 Kallang Avenue #09-10 Aperia Tower 2 Singapore 335910 PH: +65-6822-1596

DETECTION SALES, EMEA

FLIR Detection, Inc. Luxemburgstraat 2 2321 Meer Belgium PH: +32 (0) 3665 5106 detection@flir.com

www.flir.com NASDAQ: FLIR

Equipment described herein is subject to US export regulations and may require a license prior to export. Diversion contrary to US law is prohibited. Imagery for illustration purposes only. Specifications are subject to change without notice. ©2020 FLIR Systems, Inc. All rights reserved. 05/20/20

20-0749-DET-DET-Update G510 Datasheet-Ltr



The World's Sixth Sense®